


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Crystal Violet lactone - Dye content 97 %

Product Number: 332488
CAS Number: 1552-42-7
MDL: MFCD00070611
Formula: C₂₆H₂₉N₃O₂
Formula Weight: 415.53 g/mol



TEST	Specification
Appearance (Form)	Powder
Infrared Spectrum	Conforms to Structure
Wavelength	269 - 275 nm
c = 0.01g/L; Methanol	
Extinction Coefficient	≥ 50000
Carbon	72.4 - 77.8 %
Nitrogen	9.6 - 10.5 %

Remarks:

Reviewing Chemist :

Specification: PRD.0.ZQ5.10000088512

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.