

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Cyclophilin A human -  $\geq 95\%$  (SDS-PAGE), recombinant, expressed in E. coli, buffered aqueous solution

**Product Number:** C3805  
**CAS Number:** 95076-93-0  
**MDL:** MFCD00130885

**Storage Temperature:** -20 °C

TEST	Specification
mg protein/ml (UV)	$\geq 0.3$
E1% at 280nm: 4.4	
Purity by SDS Electrophoresis	$\geq 95\%$
Suitability	Pass
Tested for cleavage of suc-AAFP-PNA at 37 deg C at pH 8.0 in a coupled assay with chymotrypsin	

Specification: PRD.1.ZQ5.10000012152

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.