

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Monoclonal Anti-MDR3 P-Glycoprotein antibody produced in mouse – 250 µg/mL, clone P3 II-26, tissue culture supernatant

Product Number: M7317

MDL: MFCD03791722

Storage Temperature: -20 °C

TEST	Specification
Formulation	Conforms
Supernatant containing 0.7% BSA and 0.1% Sodium Azide	
Concentration	Conforms
250 ug/mL	
Immunocytochemistry	Conforms
1:20 - 1:50 Acetone-fixed, frozen Cytospin preparations.	
Working Dilution	Conforms
By Immunohistology; 1:20 using Acetone-fixed, frozen sections.	
Specification: PRD.0.ZQ5.10000046343	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.