


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
cis-Vaccenic acid - ≥97% (capillary GC), oil

Product Number: V0384
CAS Number: 506-17-2
MDL: MFCD00063187
Formula: C₁₈H₃₄O₂
Formula Weight: 282.46 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
Purity (GC)	≥ 97 %
% Trans Isomer	< 0.5
Impurity by GC	
Purity (TLC)	≥ 97 %

Specification: PRD.3.ZQ5.10000006183

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.