

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris(2,2,2-trifluoroethyl) borate - 97%

Product Number: **790877**

CAS Number: 659-18-7



Formula: C₆H₆BF₉O₃

Formula Weight: 307.91 g/mol

TEST	Specification
------	---------------

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Proton NMR Spectrum	Conforms to Structure
Boron NMR Spectrum	Conforms to Structure
Titration with NaOH	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000067081

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.