

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Histone Deacetylase 8 (HDAC8) antibody produced in mouse - 2.0-2.5 mg/mL, clone HDAC8-48, purified immunoglobulin, buffered aqueous solution

Product Number:

H6412

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	2.0 - 2.5 mg/ml
WB-Cell Line/Tissue Extract	Hela Cells Nuclear Fr.
WB Working Concentration	2 - 4 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000020795

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.