

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hypromellose – meets USP testing specifications

Product Number: **H3785**
CAS Number: 9004-65-3
MDL: MFCD00131360



| TEST | Specification |
|--|--------------------|
| Identity | Pass |
| Apparent Viscosity (CPS) 75% to 140% of CPS stated on the label | 3000 - 5600 |
| Loss on Drying | ≤ 5.0 % |
| Residue on ignition (Ash) | ≤ 1.5 % |
| pH | 5.0 - 8.0 |
| Methoxy Groups | 28.0 - 30.0 % |
| Hydroxypropoxy Groups | 7.0 - 12.0 % |
| Residual Solvents | Meets Requirements |

Specification: PRD.1.ZQ5.10000042764

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.