

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Zirconium(IV) silicate - -325 mesh

Product Number: 383287
CAS Number: 10101-52-7
MDL: MFCD00085353
Formula: O₄SiZr
Formula Weight: 183.31 g/mol



TEST	Specification
Appearance (Color)	Off-White
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Loss on Ignition Typically <1% Loss	
Size	≤ 44 micron

Specification: PRD.0.ZQ5.10000018990

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.