


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tropolone - 98%

Product Number: **T89702**
CAS Number: 533-75-5
MDL: MFCD00004158
Formula: C₇H₆O₂
Formula Weight: 122.12 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)

White to Yellow and Faint Beige to Beige

Appearance (Form)

Crystal to Powder

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

Specification: PRD.0.ZQ5.10000036096

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.