

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2'-Deoxyadenosine 5'-triphosphate sodium salt solution - 10 mM

Product Number: **D6920**
CAS Number: 1927-31-7

Formula Weight: 535.15 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Concentration	10 - 12 mM
DNase, Exonuclease Detection	None Detected
RNAse Detection	None Detected
NICKase, Endonuclease Detection	None Detected
Suitable for PCR	Suitable
Purity (HPLC)	≥ 99 %

Specification: PRD.1.ZQ5.10000005508

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.