

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

2-Methyltetrahydrofuran - BioRenewable, ReagentPlus<sup>®</sup>, ≥99.5%, contains 150-400 ppm BHT as stabilizer

Product Number:


155810

CAS Number:

96-47-9

MDL:

MFCD00005367



Formula Weight:

86.13 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 99.50 %

Stabilizer

150 - 400 ppm

BHT

Water (by Karl Fischer)

≤ 1.0 %

Peroxide (ppm)

≤ 100 ppm

as H<sub>2</sub>O<sub>2</sub>

Specification: PRD.3.ZQ5.1000005286

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.