


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Diethanolamine reagent grade, ≥98.0%

Product Number:**D8885**

CAS Number:

111-42-2

MDL:

MFCD00002843

Formula:

C₄H₁₁NO₂

Formula Weight:

105.14 g/mol

TEST	Specification
Appearance (Clarity)	Clear to Very Slightly Hazy
Appearance (Color)	Conforms to Requirements
Colorless to Very Faint Yellow (Viscous Liquid) or White (Solid)	
Appearance (Form)	Viscous Liquid or Solid
Infrared spectrum	Conforms to Structure
Solubility (Turbidity)	Clear
c = 50 mg/ml; Water	
Solubility (Color)	Colorless
Purity (GC)	≥98.0 %
Recommended Retest Period	-----
3 Years	

Remarks:

Specification Date : 06/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.