

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Triethylene glycol mono-11-mercaptoundecyl ether 95%

Product Number:**673110**

CAS Number:

130727-41-2

MDL:

MFCD08457662

Formula:

C₁₇H₃₆O₄S

Formula Weight:

336.53 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

Colorless to Brown

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥94.0 %

Remarks:

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.