

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:


Cholesterol - Sigma Grade, ≥99%

Product Number: C8667 CAS Number: 57-88-5

MDL: MFCD00003646

Formula: C27H46O Formula Weight: 386.65 g/mol

Storage Temperature: -20 ℃



TEST

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)

100 mg/mL, CHCl3 (APHA NMT 20)

Proton NMR spectrum Water (by Karl Fischer)

Purity (GC) Carbon

Recommended Retest Period

6 Years

Specification: PRD.2.ZQ5.10000005983

Specification

White to Off White

Pow der

Colorless to Very Faint Yellow

Clear

Conforms to Structure

< 1.0 %</pre>> 99 %

83.0 - 84.7 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.