


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
PDPP2T-TT-ODProduct Number: **791989**  
CAS Number: 1444870-74-9Formula: (C<sub>60</sub>H<sub>90</sub>N<sub>2</sub>O<sub>2</sub>S<sub>4</sub>)<sub>n</sub>

TEST	Specification
------	---------------

Appearance (Color)	Blue to Very Dark Blue and Purple to Very Dark Purple
Appearance (Form)	Powder or Solid
Infrared Spectrum	Conforms to Structure
Carbon	
Specification: Report Result	
Nitrogen	
Specification: Report Result	
Sulfur	
Specification: Report Result	
Molecular Weight by GPC	40000 - 60000
Poly Dispersity MW/MN	≤ 3.2
Product of Supplier	Confirmed
Product of Xerox(R)	

Specification: PRD.2.ZQ5.10000050077

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.