

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Bovine Serum Albumin - lyophilized powder, y-irradiated, Globulin Free, BioXtra, suitable for cell culture

 Product Number:
 A3156

 CAS Number:
 9048-46-8

 MDL:
 MFCD00130384

Formula Weight: 66,000 g/mol Storage Temperature: $2 - 8 \text{ }^{\circ}\text{C}$

TEST Specification

Appearance (Color) White to Off White

Appearance (Form) Powder

Solubility (Color) Faint Yellow to Yellow to Brown-Yellow

Solubility (Turbidity) Clear to Slightly Hazy

5g/50 mL H2O

Agarose Electrophoresis \(\geq 99 \%\)
Cell Culture Test \(\text{Pass} \)
Water (by Karl Fischer) \(< 10 \%\)

Sterilized (Gamma Irradiated) Pass

Specification: PRD.4.ZQ5.10000010746

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.