

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Dicyclopentadiene - contains BHT as stabilizer

**Product Number:** 454338  
**CAS Number:** 77-73-6  
**MDL:** MFCD00078246  
**Formula:** C<sub>10</sub>H<sub>12</sub>  
**Formula Weight:** 132.20 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless or White
Appearance (Form)	Solid or Semi-Solid or Liquid
Infrared spectrum	Conforms to Structure
GC Purity (sum of isomers) as a Mixture of Endo and Exo	≥ 96.0 %A
Stabilizer Contains BHT as Stabilizer	Confirmed

Specification: PRD.0.ZQ5.10000028763

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.