


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Terephthaloyl chloride - ≥99%, flakes

Product Number: 120871
CAS Number: 100-20-9
MDL: MFCD00000693
Formula: C₈H₄Cl₂O₂
Formula Weight: 203.02 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Flakes or Crystalline Flakes	Conforms to Requirements
Titration by AgNO ₃	99.0 - 101.0 %
GC (area %)	≥ 99.0 %
Solubility (Turbidity) Solubility Concentration: 10% in Ethanol	Clear
Solubility (Color)	Colorless

Specification: PRD.0.ZQ5.10000026699

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.