

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lithium hydride - powder, -30 mesh, 95%

Product Number:

201049

CAS Number:

7580-67-8

MDL:

MFCD00011074

Formula:

HLi

Formula Weight:

7.95 g/mol



| TEST | Specification |
|------|---------------|
|------|---------------|

Appearance (Color)

White to Grey

Appearance (Form)

Powder

Hydrogen Evolution

94.5 - 105.5 %

Size

Confirmed

-30 Mesh

Specification: PRD.0.ZQ5.10000049347

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.