

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Xylenol Orange tetrasodium salt - ACS reagent

Product Number: **398187**
 CAS Number: 3618-43-7
 MDL: MFCD00011493
 Formula: C₃₁H₂₈N₂Na₄O₁₃S
 Formula Weight: 760.58 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Very Dark Brown to Very Dark Red
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Use Test	Pass
Appearance (Clarity)	Pass
Clarity of Solution	
Meets ACS Requirements	Current ACS Specification
Recommended Retest Period	-----
3 Years	

Specification: PRD.2.ZQ5.10000025115

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.