

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Lithium tetraborate – ≥99.995% trace metals basis

Product Number:

**254304**

CAS Number:

12007-60-2

MDL:


MFCD00011083

Formula:

B<sub>4</sub>Li<sub>2</sub>O<sub>7</sub>

Formula Weight:

169.12 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Li & B Components	
Purity	Meets Requirements
> = 99.995% Based on Trace Metal Analysis	
Trace Metal Analysis	≤ 50.0 ppm

Specification: PRD.2.ZQ5.10000033856

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.