

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Hydrogen peroxide solution 30 % (w/w) in H₂O, contains stabilizer

Product Number:

H1009

CAS Number:

7722-84-1

MDL:

MFCD00011333

Formula:


H₂O₂

Formula Weight:

34.01 g/mol

Storage Temperature:

2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
% Purity (Titration)	29.0- 32.0
Recommended Retest Period 2 Years	-----

Remarks:

Specification Date : 01/25/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.