

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Trypsin from porcine pancreas - tablet, 1 mg tablet

Product Number:

T7168

CAS Number:

9002-07-7

MDL:

MFCD00082094

Storage Temperature:

-20 °C

TEST**Specification**

Total Tablet Weight	27.0 - 33.0 mg/tab
Appearance	White to Off-White
Form	Round Tablet
Diameter	5/32 inches
Dissolution Time	< 10 min
1 tablet in 2 mL H ₂ O	
pH	6.8 - 7.8
1 tablet in 2 mL H ₂ O	
Enzymatic Activity	90 - 110 %
when compared to standard	

Specification: PRD.0.ZQ5.10000046667

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.