

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

DMT-2'O-TBDMS-rC(ac)-1,3-¹⁵N₂ phosphoramidite – ≥98 atom % ¹⁵N, ≥95% (CP)

Product Number:

904562

Formula:


C₄₇H₆₄N₂N₃O₉PSi

Formula Weight:

904.08 g/mol

Storage Temperature:

-20 °C



TEST

Specification

Appearance (Form)

Powder

Purity (31P NMR)

≥ 95 %

15N Enrichment

≥ 98 atom %

Specification: PRD.0.ZQ5.10000102453

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.