


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(S,R,S)-AHPC-PEG₃-azide - ≥95%

Product Number: 904651



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow
Appearance (Form)	Liquid or Semi-Solid or Paste or Solid
¹ H NMR Spectrum	Conforms to Structure
Specific Rotation	
0.1%; DMF	
Specification: Report Only	
Carbon Content	53.0 - 58.6 %
Nitrogen Content	14.4 - 15.9 %
Purity (HPLC)	≥ 95.0 %

Specification: PRD.0.ZQ5.10000101365

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.