


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Allyl chloride

Product Number: 481378
CAS Number: 107-05-1
MDL: MFCD00000984
Formula: C₃H₅Cl
Formula Weight: 76.52 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 °C	1.4140 - 1.4170
Purity (GC)	≥ 98.0 %
Proton NMR spectra	Conforms to Structure

Specification: PRD.1.ZQ5.10000005139

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.