


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
7-Nitroindazole - ≥98%

Product Number: N7778
CAS Number: 2942-42-9
MDL: MFCD00022789
Formula: C₇H₅N₃O₂
Formula Weight: 163.13 g/mol
Storage Temperature: -20 °C



TEST

Specification

Appearance (Color)	Yellow to Brown
Appearance (Form)	Powder
Solubility (Color)	Yellow to Orange
Solubility (Turbidity)	Clear
50 mg/mL of Dimethylformamide	
Proton NMR spectrum	Conforms to Structure
Purity (TLC)	≥ 98 %
¹³ C NMR Spectrum	Conforms

Specification Date : 03/15/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.