

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lipopolysaccharides from Escherichia coli O26:B6 - purified by gel-filtration chromatography

Product Number:

L2762

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Form)	Lyophilized Powder
Appearance (Colour)	White to Faint Yellow/Tan
Solubility (Solvent)	Water
Solubility (Conc)	4.90 - 5.10 mg/ml
Solubility (Turbidity)	Faint Hazy to Hazy
Solubility (Color)	Colorless to Faint Yellow
Protein Content (Method)	Lowry-TCA
Prot. Content (% Protein)	≤ 5.00 %
Potency (Sample EU/mg)	≥ 500000 EU/mg
UV (% RNA)	≤ 20.00 %
Recommended Retest Period	-----
6 Years	

Specification: PRE.3.ZQ5.10000000870

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.