

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

trans-4-Pentylcyclohexanecarboxylic acid - 97%

Product Number:

261602

CAS Number:

38289-29-1

MDL:


MFCD00009964

Formula:

C₁₂H₂₂O₂

Formula Weight:

198.30 g/mol



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Assay

≥ 96.5 %

Specification: PRD.2.ZQ5.10000021622

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

