


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,6-Dimethylphenol - 99%

Product Number: D175005
CAS Number: 576-26-1
MDL: MFCD00002240
Formula: C₈H₁₀O
Formula Weight: 122.16 g/mol



TEST	Specification
------	---------------

Appearance (Color)

White to Yellow to Orange

Appearance (Form)

Powder or Solid or Crystals or
Crystalline Chunk(s) or Chunk(s)

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 98.5 %

Specification: PRD.0.ZQ5.10000089289

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.