

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Methylumbelliferyl sulfate potassium salt – sulfatase substrate

Product Number:

M7133

CAS Number:

15220-11-8

MDL:

MFCD00016970

Formula:


C₁₀H₇KO₆S

Formula Weight:

294.32 g/mol

Storage Temperature:

-20 °C



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Carbon

40.0 - 41.6 %

Purity (HPLC)

≥ 99 %

Impurity

≤ 0.5 %

Free 4-methyl-umbelliferone

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

5 mg/ml, H₂O

Specification: PRD.2.ZQ5.10000024828

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.