

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti- $\alpha$  Actinin antibody, Mouse monoclonal - clone BM-72.5, purified from hybridoma cell culture

Product Number:

**SAB4200813**

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgM Concentration (OD 280nm)	1.0 - 1.5 mg/ml
WB-Cell Line/Tissue Extract	Hela Cells
WB-Working Concentration	0.5 - 1.0 $\mu$ g/ml
ICC/IF Cells	HS-68 Human Foreskin Fibroblasts Cells
ICC/IF Fixation Reagent 1	Methanol
ICC/IF-Fixation Reagent 2	Acetone
ICC/IF-Working Concentration	5 - 10 $\mu$ g/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000020905

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.