


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lidocaine - powder

Product Number: **L7757**
CAS Number: 137-58-6
MDL: MFCD00026733
Formula: C₁₄H₂₂N₂O
Formula Weight: 234.34 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	70.3 - 73.2 %
Nitrogen	11.7 - 12.3 %
Purity (GC)	≥ 98 %
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity) 50 mg/ml, EtOH	Clear
Recommended Retest Period 4 Years	-----

Specification: PRD.1.ZQ5.10000049350

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.