

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Collagen from rattail Bornstein and Traub Type I (Sigma Type VII), powder

Product Number:

C8897

CAS Number:

9007-34-5

MDL:

MFCD00130825

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Solid
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear to Hazy
Few insolubles at 1 mg/mL, H ₂ O with 2 µL HOAc	
Purity by SDS Electrophoresis	Consistent with Past Lots

Remarks:

Specification Date : 02/18/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.