


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Benzoin methyl ether - 99%

Product Number: 479314
CAS Number: 3524-62-7
MDL: MFCD00008492
Formula: C₁₅H₁₄O₂
Formula Weight: 226.27 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Chunks
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
Solubility (Turbidity)	Clear
50mg/ml (5%) in Acetone	
Solubility (Color)	Colorless to Faint Yellow

Specification: PRD.0.ZQ5.10000041399

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.