

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Benzoyloxy-4-methoxybenzaldehyde - 98%

Product Number:

163953

CAS Number:

6346-05-0

MDL:


MFCD00003386

Formula:

C₁₅H₁₄O₃

Formula Weight:

242.27 g/mol

**TEST****Specification**

Appearance (Color)

Conforms to Requirements

White to Yellow to Tan

Appearance (Form)

Conforms to Requirements

Powder, Crystals or Crystalline Powder and/or Chunks

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

Specification: PRD.1.ZQ5.10000004260

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.