


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lapachol - 98%

Product Number: **142905**
CAS Number: 84-79-7
MDL: MFCD00001679
Formula: C₁₅H₁₄O₃
Formula Weight: 242.27 g/mol



TEST

Specification

Appearance (Color)	Light Yellow to Very Dark Yellow
Appearance (Form)	Powder or Crystals or Chunk(s) or Powder with Chunk(s) or Crystals with Chunk(s)
Infrared spectrum	Conforms to Structure
Purity (HPLC)	≥ 98 %
Solubility (Turbidity) 10mg/ml (1%), Ethanol	Clear
Solubility (Color)	Light Yellow to Very Dark Yellow

Specification: PRD.0.ZQ5.10000101836

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.