

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,5-Hexanediol 99% (mixture of isomers)

Product Number: **H11904**
 CAS Number: 2935-44-6
 MDL: MFCD00004557
 Formula: C₆H₁₄O₂
 Formula Weight: 118.17 g/mol



| TEST | Specification |
|----------------------------|-----------------------|
| Appearance (Color) | Colorless |
| Appearance (Form) | Viscous Liquid |
| Infrared spectrum | Conforms to Structure |
| GC Purity (sum of isomers) | ≥98.5 %A |

Remarks:

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.