

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
γ-Globulins from human blood - ≥99% (electrophoresis)

Product Number: **G4386**
CAS Number: 9007-83-4
MDL: MFCD00164089

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Very Faint Yellow to Light Yellow
Solubility (Turbidity)	Clear to Hazy
50 mg/ml in 0.9% NaCl	
Water (by Karl Fischer)	≤ 8 %
Nitrogen	14 - 16 %
% Sodium Chloride	≤ 4 %
Identity	Human Origin
Agarose Electrophoresis	≥ 99 %
Recommended Retest Period	-----
3 years	
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000025488

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.