

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Palladium(II) acetate – Vetec™ reagent grade, 98%

Product Number: **V900826**
CAS Number: 3375-31-3



Formula: C₄H₆O₄Pd
Formula Weight: 224.49 g/mol

TEST	Specification
------	---------------

Appearance (Colour)	Light Orange to Orange and Light Brown to Brown and Orange-Brown and Brown-Orange
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Palladium (Pd)	46.2 - 48.6 %

Specification: PRE.0.ZQ5.10000004046

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.