

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techsery@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methylphenidate hydrochloride

Product Number: M2892 CAS Number: 298-59-9

MDL: MFCD00058191 Formula: C14H19NO2 - HCI Formula Weight: 269.77 g/mol

TEST Specification

Appearance (Color) Appearance (Form) Solubility (Color) Solubility (Turbidity) 10 mg/mL, EtOH

Infrared spectrum

Titration Purity (HPLC) White to Off-White

Pow der Colorless Clear

Conforms to Structure 98.0 - 100.5 %

98.0 - 102.0 %

Specification: PRD.0.ZQ5.10000021840

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.