

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium fluoride – anhydrous, powder, 99.99% trace metals basis

Product Number: 450022
CAS Number: 7681-49-4
MDL: MFCD00003524
Formula: FNa
Formula Weight: 41.99 g/mol

NaF

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Gravimetric Analysis % F With Lead Acetate	
ICP Major Analysis Confirms Sodium Component	Confirmed
Purity 99.99% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 150.0 ppm

Specification: PRD.0.ZQ5.10000042740

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.