

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Neodymium - ingot (in mineral oil), 99.9% trace rare earth metals basis

Product Number:**261157**

CAS Number:

7440-00-8

MDL:

MFCD00011130

Formula:

Nd

Formula Weight:

144.24 g/mol

Nd

TEST	Specification
Appearance (Color) Grey	Conforms
Appearance (Form) Ingots in Oil	Conforms
ICP Major Analysis Confirms Neodymium Component	Confirmed
Purity 99.9% Based On Rare Earth Analysis	Meets Requirements
Trace Rare Earth Analysis	≤ 1500.0 ppm

Specification: PRD.2.ZQ5.10000005385

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.