


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium hexachloropalladate(IV) - 99%

Product Number: 334502
CAS Number: 16919-73-6
MDL: MFCD00043081
Formula: Cl₆K₂Pd
Formula Weight: 397.33 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Red-Orange to Red	
Appearance (Form)	Powder or Crystals
Gravimetric Analysis	26.3 - 27.2 %
% Pd	
X-Ray Diffraction	Conforms to Structure

Specification: PRD.0.ZQ5.10000033435

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.