

Product Specification

Product Name:

LL-37 (human) trifluoroacetate salt - $\geq 95\%$ (HPLC)

Product Number:

94261

CAS Number:

597562-32-8

Formula Weight:

4.492,28 g/mol

Storage Temperature:

-20 °C

| TEST | Specification |
|----------------------------|--------------------|
| Appearance (Form) | Powder |
| Appearance (Color) | Colorless or White |
| Solubility (Solvent) | Water |
| Solubility (Conc) | 1.00 - 1.04 mg/ml |
| Solubility (Color) | Colorless |
| Solubility (Turbidity) | Clear |
| Assay (Elemental Analysis) | $\geq 70.0\%$ |
| Purity (HPLC Area %) | $\geq 95.0\%$ |

Specification: PRE.1.ZQ5.10000014258

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

