

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans-2-Nonen-1-ol - ≥96%, FCC, stabilized, FG

Product Number: W337900
CAS Number: 31502-14-4
MDL: MFCD00014058
Formula: C₉H₁₈O
Formula Weight: 142.24 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.444 - 1.452
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 96.0 %
Specific Gravity	0.830 - 0.850
At 25 Deg C	
Alpha-Tocopherol	0.01 - 0.20 %
Solubility (Turbidity)	Clear to Very Slightly Hazy
c = 1ml/1ml; 95% Ethanol	
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.3.ZQ5.10000008134

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.