

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Casein Blocking Buffer 10x - for Northern and Southern blotting, solution

Product Number:

B6429

TEST	Specification
Appearance (Turbidity)	Hazy
Appearance (Color)	White to Off-White
Appearance (Form)	Liquid
Suitability	Suitable
Suitable as a blocking buffer to decrease the signal generated by non-specific binding.	

Specification: PRD.0.ZQ5.1000070176

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.