

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,5-Pentanedithiol - 96%

Product Number: 242551
CAS Number: 928-98-3
MDL: MFCD00004908
Formula: C₅H₁₂S₂
Formula Weight: 136.28 g/mol

HSCH₂(CH₂)₃CH₂SH

TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 95.5 %

Specification: PRD.1.ZQ5.10000050780

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.