

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Malt Extract Agar – Microbiologically tested

Product Number: **M6907**

TEST	Specification
Appearance (Color)	Off White to Yellow to Beige
Appearance (Form)	Powder
Solubility (Color)	Colorless to Yellow to Brown-Yellow
Solubility (Turbidity)	Clear to Hazy
3.36 g/100 mL, H ₂ O	
pH	4.5 - 4.9
(3.36% Solution)	
Performance	Pass
Expiration Dated Product	-----

Specification: PRD.5.ZQ5.10000016884

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.