

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Freund's Adjuvant, Complete - cell suspension

Product Number: **F5881**

MDL: MFCD00131105

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Light Yellow to Yellow
--------------------	------------------------

Appearance (Form) with particulates	Liquid
--	--------

Appearance (Turbidity)	Clear to Slightly Hazy
------------------------	------------------------

Emulsification	Pass
----------------	------

Forms emulsion with 0.85% NaCl

Note

Each mL contains 1 mg Mycobacterium tuberculosis(H 37RA, ATCC 25177), heat killed and dried, 0.85 mL paraffin oil and 0.15 mL mannide monooleate.

Specification: PRD.1.ZQ5.10000016267

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.